

Kapsafe™ Recapper FAQs

- 1- What is KapSafe's foot print?

60 wide x 60 deep x 83cms tall.
- 2- What is KapSafe's throughput?

Up to 20 tubes per minute, making KapSafe valuable for the large and small volume labs.
- 3- What rack types can be accommodated?

Abbott Architect, Beckman AU, Roche Cobas, Siemens Centaur XP, LGP 10 position racks.

All racks listed, except Architect racks, may be loaded intermixed on KapSafe.
- 4- What tube sizes can be run simultaneously?

All standard original vacuum collection tubes, glass or plastic, from 13-16mm in diameter 75-100mm tall, as well as 12x75mm transfer tubes. See picture of tubes on page 2.
- 5- Can KapSafe skip empty positions within a rack?

Absolutely, empty positions within a rack will be skipped to the next available tube.
- 6- What is the cap hopper capacity?

3,600 caps - KapSafe uses a **unique** cap. See picture on page 2.
- 7- Does KapSafe have continuous load and unload capability?

Yes, both the input and output bays allow for continuous loading and unloading thus maximizing the efficiency and productivity of KapSafe.
- 8- Is there a warranty and service contracts?


Yes, there is a one year warranty and service contracts available after the warranty
- 9- Is a compressor needed?

No compressed air is needed.
- 10- What's the power consumption?

Fused at 1 Amp and rated at 36 watts, KapSafe plugs into any regular outlet.
Available in 110 VAC 60 Hz & 220 VAC 50 Hz



Pictures of the KapSafe caps and tubes.



11- What are the LGP generic racks?

These are 10 position racks that accommodate 13x75mm and/or 13x100mm tubes.

12-

The LGP generic racks can be stored away in your refrigerator for archiving as they are low cost, compact and user friendly.....



PLACE BAR CODE LABELS on racks/trays for easy retrieval of samples

Furthermore, when it is time to discard the sample tubes, just flip the LGP tray and all 100 samples are disposed of into the bio-hazardous waste bin at once.

