

# CLINICON<sup>®</sup>

Laboratory Consumables & Equipment

## Introducing The New KapSafe™ Mini Benchtop Recapper

Safely recap your sample tubes with a high-speed pneumatics free benchtop recapper. Recaps original vacuum tubes 13-16mm diameter, 75-115mm tall. KapSafe accommodates most analyser racks and provides single rack recapping.



### Technical data

Capacity:	Recaps up to 1,000 tubes per hour
Sample tubes:	Standard vacuum blood tubes
Caps:	Uses cap 2174-3056
Size:	35 x 55 x 66cms (Depth x Width x Height)
Voltage:	230VAC 50Hz / 110VAC 60Hz
Power:	60VA