

Pluggo Architect 5 Position Automated Decapper - Technical Description

Pluggo is designed for the safe removal of original caps from most standard vacuum collection tubes. Pluggo can handle 5ml, 7ml and 10ml vacuum tubes with 75 x 13 mm, 100 x 13 mm, 100 x 16 mm tubes carried in the Architect racks. A minor modification to the racks is necessary when using 5ml tubes.

PLUGGO



A detachable chute guides the caps into a disposal area after decapping.

Technical data

Capacity:	Up to 5 samples in each rack
Sample tubes:	Standard blood collection tubes
Speed:	2 seconds per rack
Size:	29 x 27 x 21cm (Depth x Width x Height)
Weight:	8,0 kg
Voltage:	230VAC 50Hz / 110VAC 60Hz
Power:	60VA