

Pluggo Beckman AU 10 Position Automated Decapper - Technical Description

Pluggo is designed for the safe removal of original caps from most standard vacuum collection tubes. Pluggo can handle 5ml, 7ml and 10ml vacuum tubes with 75 x 13 mm, 100 x 13 mm, 100 x 16 mm tubes carried in the Beckman 10 position rack.

PLUGGO



A detachable chute guides the caps into a disposal area after decapping.

Technical data

Capacity:	Up to 10 samples in each rack
Sample tubes:	Standard vacuum blood tubes
Speed:	20 seconds per rack
Size:	42 x 38 x 21cm (Depth x Width x Height)
Weight:	9,0 kg
Voltage:	230VAC 50Hz / 120VAC 60Hz
Power:	60VA