

Pluggo OCD (Vitros) Automated Decapper – Technical description

Pluggo is designed for safe removal of original caps from BD Vacutainer™ vacuum collection tubes. Pluggo handles Vacutainer™ vacuum tubes with 13 x 75 mm and 13 x 100 mm size tubes carried in OCD universal 10-position carousels.

PLUGGO



The caps are guided to an outlet for collection of biohazard disposal with a detachable chute.

Technical data

Capacity:	Up to 10 samples in each rack
Sample tubes:	BD Vacutainer™
Speed:	Approx. one second per tube
Size:	430x270x210 mm (Depth x Width x Height)
Weight:	8,0 kg
Voltage:	230VAC 50Hz / 110VAC 60Hz
Power:	60VA